

THE "2 IN 1" HOOF

The HOOF-it® 'All in One' Hoof Stand, designed by Farrier Steve Samet, features a shock absorbing, integrated rubber cradle and hoof post into one easy to use unit.

The Hoof Stand provides stable and safe support for the horse, allowing it to relax without putting all of its weight on you. For additional safety and stability when working with the Hoof-it® Stand try the following:

- * keep one foot on the base unit when working with the Hoof-it® Hoof Stand
- * do not leave the horses' leg or hoof in the cradle or on the stand when unattended.

Benefits for Farriers, Veterinarians and Horse Owners who use the HOOF-it® Hoof Stand include:

- √ Easy to use "All in One" innovative unit design integrates both the post and cradle into one.
- √ The innovative design of the HOOF-it® Hoof Stand dramatically reduces your risk of getting injured during hoof care tasks.
- √ The HOOF-it Stand is ideally suited for use with all types of horses, including older horses that require special support and comfort during hoof care tasks.
- √ The rubber cushioned cradle fits different hoof sizes and hugs the hoof.

Let HOOF-it® Hoof Stand deal with the majority of the horses weight during any type of hoof care tasks, rather than your back or knees.




From Stand to Cradle in one easy step.



- Working heights (measured from the ground):
- Stand:
 - low: Horse: 20" / 51 cm - Draft: same
 - fully extended: Horse: 27" / 68 cm - Draft: same
 - Standard cradle:
 - top width: Horse: 6 1/2" / 16.5 cm - Draft: 9 1/2" / 24 cm
 - depth: Horse: 2 1/2" / 5.5 cm - Draft: 3 1/4" / 7.8 cm
 - Base:
 - High (ion top of the clip): Horse: 13" / 33 cm - Draft: same
 - Width : Horse: 18" / 46 cm - Draft: same
 - Ball :
 - height: Horse: 1 5/8" / 4.5 cm - Draft: same
 - diameter 2 1/2" / 6.5 cm - Dract: same

ASSEMBLE

1. A) Stand B) Cradle and C) Green Base.



2. Take A (stand) and insert into B (cradle)



3. Put both A (stand) and B (Cradle) into C (green base).



4. After A (stand) and B (cradle) are inserted into C (green base) you will need to make a one time tension adjustment. Tighten the clamping mechanism tight enough so that when clamp is in the "locked" position it will hold part A (the stand) up without any slippage. Do this by pulling the clamp handle out straight (in a loose or "unlocked" position) then, while holding the left nut, turn the handle clockwise until you hit the desired tension (test by pulling the hoof stand (part A) up and lock clamp, making sure it is holding the stand up tightly). When you have found your desired tension, it will be one easy step to switch from stand to cradle and adjust the height of your stand or cradle.



The HOOF-it Hoof Stand is the only stand on the market that switches from a stand to a cradle in one easy step. NO need to remove any parts. The stand retracts right into the cradle.

Made for the professional farrier, equine owner and Veterinarian.

HOOF-it® Technologies

PO Box 382

Minden, Nevada 89423

Toll Free: (888)665.9099

Phone: (775)783.9099

Fax: (775)783.1011

www.hoof-it.com



HOOF-it "2 in 1"
Hoof Stand



www.hoof-it.com